

Ready for life transfer aids

The Ready for life transfer aids are clever and easy-to-use assistive devices for safe and easy transfers of patients. They offer cost effective solutions to many challenging transfer situations. At the same time they spare the backs and shoulders of care personnel and prevent work related injuries which, in turn, achieves significant cost savings. The product category include RoMedic's classic manual Ready-products as well as a completely new mobile and collapsible patient lift. The transfer assistive products can be used either independently or as a complement to products in the SystemRoMedic range.

ReadyGlide

ReadyGlide is a sliding board that has been designed to simplify transfers in a seated position, for example, between bed and wheelchair or between wheelchair and toilet or car.



Item No.	Product	Measurements
15030	ReadyGlide, polyethylene	60 x 25 cm



ReadySlide

ReadySlide is a thin, pliable and very easy-to-use sliding mat that facilitates various types of patient transfers.




Item No.	Product	Measurements
11250	ReadySlide, polyamid	70 x 70 cm
11230	ReadySlide, polyamid	110 x 70 cm
11235	ReadySlide, polyamid, extra wide	110 x 90 cm
11210	ReadySlide, polyamid	195 x 70 cm
11215	ReadySlide, polyamid, with handles	195 x 70 cm



ReadyLett

ReadyLett is a mattress cover that is used in situations where there are special requirements for patient mobility in bed, for example, for bedridden patients in intensive-care units.




Item No.	Product	Measurements
11480	ReadyLett, polyamid	200 x 90 cm



ReadySheet

ReadySheet is a simple and versatile assistive device for repositioning in one location, for example, for turning or moving the user higher up in bed.




Item No.	Product	Measurements
11440	ReadySheet, polyamid	175 x 75 cm
11445	ReadySheet, polyamid, tubular	175 x 75 cm
11412	ReadySheet, polyamid	180 x 110 cm
11415	ReadySheet, polyamid, package with 2 pcs 11412	180 x 110 cm
11450	ReadySheet, polyamid	150 x 75 cm
11455	ReadySheet, polyamid, package with 2 pcs 11450	150 x 75 cm
11430	ReadySheet, polyamid	195 x 75 cm
11425	ReadySheet, polyamid	200 x 90 cm



ReadyMove

With ReadyMove, using minimal effort and without lifting, a caregiver can give a user safe and secure transfer and mobility support.



Item No.	Product	Measurements
13010	ReadyMove, nylon	20 x 50 cm



ReadyBelt

ReadyBelt is a belt that is designed to provide the user with support and help during mobility training and different transfer situations. For the caregiver, ReadyBelt offers an ideal means of supporting the user during mobility training or visits to the toilet.



Item No.	Product	Measurements
16023	ReadyBelt, polyester. Small, 7 handles	105 x 10 cm
16024	ReadyBelt, polyester. Medium, 9 handles	120 x 10 cm
16025	ReadyBelt, polyester. Large, 11 handles	140 x 10 cm



ReadyGrip

ReadyGrip is a very easy-to-use assistive device that facilitates raising from a recumbent position to a seated position. ReadyGrip is intended for active patients who have arm strength and good torso stability.



Item No.	Product	Measurements
16090	ReadyGrip, nylon	170 x 17 cm



Carina350EM/EE

A versatile and priceworthy solution for safe and easy lifting

Carina350 is a collapsible, mobile patient lift adapted for safe and easy lifting of users weighing up to 165 kg/350 lbs. Carina350 is available in two models, with manual and electrical base widening, respectively. Both models are sturdy, stable and safe and, at the same time, light-weight for steel lifts. Quite simply, Carina350 offers a unique combination of functionality, performance and quality in relation to price. Carina350 is delivered complete and ready to use and is equipped with a stable steel sling bar with reliable safety latches.



Item No.	Product
60600011	Carina350EM, mobile lift with manual base widening
60600001	Carina350EE, mobile lift with electrical base widening





All you need for safe and easy lifting of users weighing up to 230 kg/505 lbs

Reda is a safe and easy-to-use ceiling lift unit that combines new design with well-thought-out and proven construction and functionality. Reda can be used on all types of rails, both permanent rail systems and on portable, freestanding lift stands. Reda also functions equally well in all settings and in all types of rooms. It is very easy to operate, moves silently and gently and features many functions that maximize comfort and security for both the user and the caregiver. Reda is delivered complete with a sling bar and a new, uniquely user-friendly hand control.


Item No.	Product
50100023	Reda, ceiling lift unit with sling bar and hand control



Milky Way

A simple and safe rail system that makes the most of your specific conditions

Milky Way is a complete, ingeniously simple and safe rail system for permanent installation on the ceiling and/or on the walls. Each rail system is individually designed in close cooperation with the customer, allowing adaptation for specific conditions and to meet specific requirements. Together with an optional choice of lift unit and the right lifting accessories, Milky Way offers unlimited possibilities for developing and adapting the perfect stationary lift solution for all settings, rooms and lifting needs. For more information about the Milky Way rail system, please contact RoMedic AB or your local Ready for life distributor. With your conditions and requirements as a starting point, we will help you create the perfect lifting solution for you.



Ready for life

Ready for life – a new product series for improved quality of life for all

No one life is like the other. We all have different living conditions and during a lifetime, these conditions are constantly changing. Good and bad things happen to us all and, if we are lucky, we all get older, day by day.

Life, whatever it looks like, is a gift, not a given. And it should not just be about existing or surviving. Life should be lived. We don't have to be young, healthy or even independent to live. But we do have to be ready, that is willing and able, for life.

Functional disabilities, whether they are permanent or temporary conditions, ideally should not get in the way of anyone's ability to live a full and satisfactory life. The will comes from within, and when it comes to ability, Ready for life can improve the conditions by providing tools to better enable all people to live a fuller, more active life, every day of their lives.

Ready for life is a series of assistive devices that are designed to give users the best possible conditions for living an active life on a daily basis. The products have been specially developed to meet the conditions for home-care and institutional care settings. Ready for life offers the right quality, the right function and the right performance at the right price.

Ready for life transfer aids

The Ready for life transfer aids are clever and easy-to-use assistive devices for safe and easy transfers of patients. They offer cost effective solutions to many challenging transfer situations. At the same time they spare the backs and shoulders of care personnel and prevent work related injuries which, in turn, achieves significant cost savings. The product category include RoMedic's classic manual Ready-products as well as a completely new mobile and collapsible patient lift. The transfer assistive products can be used either independently or as a complement to products in the SystemRoMedic range.



© RoMedic AB - 100517 - www.romedic.com

RoMedic

A company in the Handicare Group

RoMedic AB

Box 640 • Veddestavägen 15, SE-175 27, Järfälla, Sweden
 Tel: +46 8 594 773 33 • Fax: +46 8 594 773 44
 www.romedic.com • info@romedic.com