

Accessories for lifting slings



AdjustmentLoop™ SystemRoMedic

AdjustmentLoop extension straps are used to extend the sling straps. They are available in two lengths and are delivered in pairs.

Item No.	Product	Description
49182006	AdjustmentLoop	Straps 10 cm/4", pair
49182007	AdjustmentLoop	Straps 20 cm/8", pair



SoftLegSupport™ SystemRoMedic

Pressure-sensitive users sometimes need extra soft leg supports. SoftLegSupport is slipped on each leg support of the sling. Delivered in pairs.

Item No.	Product	Description
49373003	SoftLegSupport	Padding for leg supports XS, pair
49373006	SoftLegSupport	Padding for leg supports S/M, pair
49373007	SoftLegSupport	Padding for leg supports L/XL, pair



StandUpSling™ Padding SystemRoMedic

StandUpSling can be equipped with extra padding as an accessory for pressure-sensitive users.

Item No.	Product	Size
45273306	StandUpSling Padding	Fits StandUpSling size S-M
45273307	StandUpSling Padding	Fits StandUpSling size L-XL



HeadSupport™ SystemRoMedic

HeadSupport is used as a complement for slings without head support. HeadSupport is placed around the head and attached to the sling bar.

Item No.	Product	Description
49302000	HeadSupport	For extra head support



NameTag™ SystemRoMedic


NameTag is a textile marking label which is easily attached to, removed from and moved between products. It is safely attached and can be washed at up to 90°C/194°F.

Item No.	Product	Description
49200006	LimbSling	Lift strap in polyester
49292000	LimbSlingDisposable	Lift strap in disposable material, 10-pack
49700000	NameTag	Marking label, 10-pack

Item No.	Product	Description
49200006	LimbSling	Lift strap in polyester
49292000	LimbSlingDisposable	Lift strap in disposable material, 10-pack

* LimbSlingDisposable has no pocket on the strap

Other lifting accessories



StretcherBar™ SystemRoMedic

StretcherBar is a slingbar used together with StretcherSling for horizontal lifting. It can be used with, among others, the Rise ceiling lifts and Eva mobile lifts. The width of the sling is easily adjusted by pulling out or pushing in the suspension hooks. The equilibrium is easily set with a built-in control so that the user is lifted horizontally. The built-in control means that no lifting height is lost. The slingbar is made from aluminium, giving many environmental advantages.

Item No.	Product
70200006	StretcherBar 10 point
46502007	StretcherSling

Safe working load: 300 kg/660 lbs


Both Scale660 and Scale880 are used together with Rise ceiling lifts and Eva mobile lifts for weighing of users.

Clear display.

Item No.	Product
70100002	Scale660
70100003	Scale880

Safe working load: 400 kg/880 lbs

Scale880™ SystemRoMedic



AmbulationArm™ SystemRoMedic

The armrests can be employed by the user during gait training but they can also be used by the carer for manoeuvring the lift.

Item No.	Product	Description
80100010	AmbulationArm	Armrests - fits the Eva lifts

SystemRoMedic™ "We love easy transfers"

The philosophy behind SystemRoMedic is simple. We love to develop and manufacture assistive devices that make life easier for both users and personnel in the care sector.

SystemRoMedic comprises four product areas that meet most needs in different transfer situations.

- **Transfer**, between two locations.
- **Positioning**, within the same location.
- **Support**, for mobility.
- **Lifting**, both manual and mechanical.



SlingBar™ SystemRoMedic

Slingbars with three different widths: 350, 450, 600 mm/13.7", 17.7", 23.6"

SlingBar in aluminium is available in different widths for users weighing up to 300 kg/660 lbs. The sling bars have safety latches which prevent the sling straps from creeping out of the sling bar.

SystemRoMedic's standard SlingBar is 450 mm/17.7" wide and excellent for most lifting situations. It is also available with most of SystemRoMedic's lifts.

For safety reasons, it is important to choose the right sling bar for the right sling and user. A narrow sling bar suits children and small adults for normal lifting but it is also suitable for gait training with a WalkingVest. A wide sling bar is excellent for a large user with increased need for room within the sling.

To get more room in a two-point sling bar, you can use SideBar sidebars as an accessory. They open up the sling and provide for a more reclined position.

Read the instruction manuals for the slings for recommendations on which sling bar is best suited.



SideBar, fits Rino600.

Another accessory for SlingBar is the SwiftHook. With SwiftHook changing sling bars or connecting a scale becomes very easy.



SwiftHook.

SlingBar can also be equipped with a soft protection, SlingBarWrap, which is sometimes needed for users with uncontrollable movements. It is attached with Velcro and can easily be removed for cleaning. There is also a SlingBarWrap for MiniLift200.



SlingBarWrap.

Item No.	Product	Description
70200001	SlingBar S	350 mm/13.7", two-point slingbar
70200002	SlingBar M	450 mm/17.7", two-point slingbar
70200003	SlingBar L	600 mm/23.6", two-point slingbar
70200013	SideBar	Sidebars for Rino600, pair
70200008	SwiftHook	Quick-release hook
70200009	SlingBarWrap	Soft protection for slingbar S, M, L
70200010	Room Transfer Hook (RT-H)	Quick connection "room to room"
70200011	H-adapter	Adapter for 70200010
70200012	SlingBarWrap MiniLift	Soft protection for MiniLift 200

RoMedic™ A company in the Handicare group

RoMedic AB
Veddestavägen 15, Box 640, 175 27 Järfälla, Sweden
Tel. +46 (0) 8 594 773 33 • Fax. +46 (0) 8 594 773 44
www.romedic.com • info@romedic.com

Lifting slings and accessories SystemRoMedic



The lifting sling makes the difference

SystemRoMedic offers a wide variety of functional, high quality lifting slings adapted to meet the different requirements of all kinds of lifting situations and users. The lifting slings are available in four different materials and in sizes ranging from XXS to XXL. All lifting slings can, to an advantage, be used together with SystemRoMedic's mobile lifts and ceiling lift systems. Naturally they can also be used together with lifts of other makes.

We love easy transfers and, when developing lifting slings and other accessories for our lifts, we are always looking for the very best materials and the simplest solutions for both carer and user. Our close cooperation with carers and users is highly contributing to the fact that our products are both safe and very easy to use and that they, at the same time, provide the user with both security and comfort. All SystemRoMedic lifting slings have a safe working load of up to 300 kg/660 lbs.

Consider this when selecting sling size and material.

Size

It is too small

- when the sling bar gets too close to the user's face
- when the leg supports cut into the user's crotch
- when the sling looks too tight for the user

It is too large

- when the user's buttocks slide out too much through the opening
- when the leg supports slide forwards towards the back of the user's knees

Material

Polyester

- Wear resistant
- Low friction
- Withstands high washing temperature

Polyester net

- Airy
- Water runs through
- Withstands high washing temperature
- Can be left under the user

Plastic-coated net


- Does not absorb water, good for bathing
- Can be disinfected with alcohol

Disposable material

- Strong disposable material
- Is not cleaned but disposed of after use
- Prevents the spread of multi-resistant bacteria

Important details

SystemRoMedic has a comprehensive assortment of well made, comfortable and functional slings with many practical details. Always read the manual before using the sling and don't forget to perform regular function tests. Read more on www.romedic.com



Lifting slings

Lifting between bed and wheelchair



LowBackSling
SystemRoMedic

LowBackSling is a sling which gives the user support from under the legs to up under the arms. LowBackSling is designed for use in the most common lifting situations and when the user has trunk stability.

Product specification	Item No.
Polyester XS-S	4020003-4
Polyester M-XL	4020006-8



BasicSling
SystemRoMedic

BasicSling is a sling which gives the user support from under the legs to up to the shoulders. BasicSling is designed for use in the most common lifting situations and when the user has reduced trunk stability but good head stability.

Product specification	Item No.
Polyester XS-S	4050003-4
Polyester M-XL	4050006-8
Polyester net XS-S	4052003-4
Polyester net M-XL	4052006-8



HighBackSling
SystemRoMedic

HighBackSling is a sling which is easy to fit for the carer, while providing security for the user. HighBackSling is designed to meet a need for considerable support.

Product specification	Item No.
Polyester XS-S	4080003-4
Polyester M-XXL	4080006-9
Polyester net XS-S	4082003-4
Polyester net M-XXL	4082006-9
Plastic-coated net XS-S	40852003-4
Plastic-coated net M-XXL	40852006-9



ShellHBSling
SystemRoMedic

ShellHBSling is a sling which is easy to fit and provides security for the user. ShellHBSling is slightly cup-shaped after the body and thereby provides for an upright seated posture and excellent support for the head.

Product specification	Item No.
Polyester XS-S	4070002-4
Polyester M-XL	4070006-8

Lifting from floor



HighBackSling
SystemRoMedic

HighBackSling is a sling which is easy to fit and provides security for the user. HighBackSling is designed to provide considerable support and is excellent for lifting from the floor.

Product specification	Item No.
Polyester XS-S	4080003-4
Polyester M-XXL	4080006-9
Polyester net XS-S	4082003-4
Polyester net M-XXL	4082006-9
Plastic-coated net XS-S	40852003-4
Plastic-coated net M-XXL	40852006-9



ShellHBSling
SystemRoMedic

ShellHBSling is a sling which is easy to fit and provides security for the user. ShellHBSling is slightly cup-shaped after the body and thereby provides for an upright seated posture and excellent support for the head.

Product specification	Item No.
Polyester XS-S	4070002-4
Polyester M-XL	4070006-8



FlexibleSling
SystemRoMedic

FlexibleSling is a sling providing security while simultaneously being soft and flexible around the body. FlexibleSling fulfills the need for a lot of support while, at the same time, adapting naturally to the contours of the body. It can be left in a wheelchair, moulded seat etc.

Product specification	Item No.
Polyester net XS-S	42522002-4
Polyester net M-XL	42522006-8
Polyester net (m.first) XS-S	42520002-4
Polyester net (m.first) M-XL	42520006-8

To/from toilet



HygieneSling
SystemRoMedic

HygieneSling is a sling which facilitates dressing/undressing of the user in connection with toileting. HygieneSling requires the user to have a sufficient degree of tonus.

Product specification	Item No.
Polyester (with belt) XS-S	43504003-4
Polyester (with belt) M-XL	43504006-8
Polyester S	43604004
Polyester M-L	43504006-7



HygieneHBSling
SystemRoMedic

HygieneHBSling is a sling which gives a lot of support. The seated posture is upright and the opening around the buttocks is adapted to enable removal of underwear while the user is sitting in the sling.

Product specification	Item No.
Polyester XS-S	43704003-4
Polyester M-XL	43704006-8



StandUpSling
SystemRoMedic

StandUpSling is a vest designed for use in combination with a sit-to-stand lift. It is easy to fit for the carer, provides the user with support and an extra feeling of security and also facilitates toileting.

Product specification	Item No.
Polyester S	45200004
Polyester M-XL	45200006-8

Accessories:

Product specification	Item No.
StandUpSling Padding S-M	45273306
StandUpSling Padding L-XL	45273307



ThoraxSling
SystemRoMedic

ThoraxSling is a comfortable sit-to-stand vest used together with MiniLift200 and MiniLift125. It is easy to fit for the carer, provides the user with support and a feeling of security and facilitates toileting.

Product specification	Item No.
ThoraxSling S	691591
ThoraxSling M-XXL	691592-5
ThoraxSling (w. SeatSupport) S	691596
M-XL	691597-9
CallStrap	691174
Attachments for CallStrap 1 pair	400962000

*Attachments are necessary for use

Gait training



WalkingVest
SystemRoMedic

WalkingVest is a vest designed for gait training with the help of a stationary or mobile lift. WalkingVest provides support around the trunk and is kept in place with a padded leg harness.

Product specification	Item No.
Polyester XS-S	4420002-4
Polyester M-XXL	4420006-9

Rising/standing



StandUpSling
SystemRoMedic

StandUpSling is a sit-to-stand vest designed for use with a lift. It is easy to fit for the carer and provides the user with an extra feeling of security. StandUpSling also facilitates toileting.

Product specification	Item No.
Polyester S	45200004
Polyester M-XL	45200006-8

Accessories:

Product specification	Item No.
StandUpSling Padding S-M	45273306
StandUpSling Padding L-XL	45273307



ThoraxSling
SystemRoMedic

ThoraxSling is a comfortable sit-to-stand vest used together with MiniLift 200 and MiniLift125. It is easy to fit for the carer, provides the user with support and a feeling of security and facilitates toileting.

Product specification	Item No.
ThoraxSling S	691591
ThoraxSling M-XXL	691592-5
ThoraxSling (w. SeatSupport) S	691596
M-XL	691597-9
CallStrap	691174
Attachments for CallStrap 1 pair	400962000




StandingVest
SystemRoMedic

StandingVest is a lift vest used for assisting a user using a standing shell. The vest provides support in the crotch, around the trunk and all the way up behind the head.

Product specification	Item No.
Polyester XS-S	44503003-4
Polyester M-XL	44503006-8

Bath



HighBackSling
SystemRoMedic

HighBackSling bath net is excellent for bathing. HighBackSling is a sling which is easy to fit for the carer, while providing security for the user. HighBackSling is designed to meet a need for considerable support. Another alternative for bathing can be a sling in polyester net.

Product specification	Item No.
Plastic-coated net XS-S	40852003-4
Plastic-coated net M-XXL	40852006-9



StretcherSling
SystemRoMedic

StretcherSling is intended for horizontal lifting together with StretcherBar. It is easy to fit for the carer, while simultaneously providing security for the user.

Product specification	Item No.
StretcherSling Polyester	46502007
StretcherBar	70200006



SafeHandlingSheet
SystemRoMedic

SafeHandlingSheet is used for lifting and positioning of a user in the bed without loading the carer. SafeHandlingSheet is used together with a ceiling or mobile lift. The loops of the SafeHandlingSheet are attached to a regular sling bar.

SafeHandlingSheet is available in a disposable material which is disposed of after use instead of being cleaned. This prevents the spread of multi-resistant and other types of contagious bacteria and viruses.

Product specification	Item No.
SafeHandlingSheet, fourteen lift straps (seven on each side) with three alternative lifting loops	47592006



SafeHandlingSheet is also available in polyester net. This sheet can be left in the bed, since it is a strong and airy material to lie on.

Product specification	Item No.
SafeHandlingSheet, sixteen lift straps (eight on each side) with three alternative lifting loops	47522007

Disposables



DisposableSling
SystemRoMedic

DisposableSling is made from a disposable material which is disposed of after use instead of being cleaned. This prevents the spread of multi-resistant and other types of contagious bacteria and viruses. DisposableSling is designed for users needing an upright sitting position with support for the legs and behind the head. It has leg supports for easy application and removal. The user has the arms on the inside. Made in several sizes.

Product specification	Item No.
HighBack S, 10-pack	40992004
HighBack M, 10-pack	40992006
HighBack L, 10-pack	40992007
HighBack XL, 10-pack	40992008



SafeHandlingSheet
SystemRoMedic

SafeHandlingSheet is available in a disposable material which is disposed of after use instead of being cleaned. This prevents the spread of multi-resistant and other types of contagious bacteria and virus. SafeHandlingSheet is used for lifting and positioning users in the bed without loading the caregiver.

LimbSling
SystemRoMedic

LimbSling is used for lifting single body parts for treatment and wound cleaning. LimbSling is available in a strong disposable material and disposed of after use instead of being cleaned.

Product specification	Item No.
Disposable 10-pack, 300 kg/660 lbs	49292000

